


Molecular carbapenemase tests Eazyplex SuperBug Complete B, Novodiag CarbaR+, and Amplidiag CarbaR+MCR, provide high sensitivity in the detection of carbapenemase-positive bacteria



Tanja Holma, Jenni Antikainen, Anu Pätäri-Sampo, and Johanna Haiko

Contact: johanna.haiko@hus.fi

Helsinki University and Helsinki University Hospital, HUSLAB, Department of Clinical Microbiology, Helsinki, Finland

Introduction

- Carbapenemase-producing *Enterobacteriaceae* (CPE), MDR *Pseudomonas aeruginosa*, and MDR *Acinetobacter*: rapid screening needed
- The new regulations by the European Union regarding in vitro diagnostic medical devices → replacement of in-house methods with commercial tests
- Aim: to compare three CE-IVD marked carbapenemase tests:
 - Eazyplex SuperBug Complete B (Amplex Diagnostics)
 - Novodiag CarbaR+ (Mobidiag)
 - Amplidiag CarbaR+MCR (Mobidiag)

Materials and methods

- Workflow and the basic characteristics of the tests: see Fig. 1
- Material (Table 1):
 - Clinical bacterial isolates (n=32) and quality control strains (n=5)
 - Clinical samples (n=34), mostly from rectum, urine, or wounds
 - Some strains/samples included several target genes
- Methods:
 - One colony
 - in RALF for Eazyplex
 - in eNAT (Copan) for Novodiag and Amplidiag
 - Clinical swabs or 300 µl of liquid samples in eSwab (Copan)
 - as such for Eazyplex
 - 200 µl into an eNAT for Novodiag and Amplidiag
 - Comparison to traditional culture:
 - mSuperCARBA (CHROMagar) plate overnight
 - ID with MALDI-TOF MS (bioMérieux)
 - in-house carbapenemase PCR (Pasanen et al., 2014)

Figure 1. Workflow chart




Table 1. Bacterial strains and samples

Target genes	Clinical bacterial isolates (n=32)	Quality control strains (n=5)	Clinical samples (n=34)
KPC	4	0	6
NDM	3	1	7
IMP	2	1	0
VIM	2	1	0
OXA-48/181	6	2	10
MCR	0	1	0
GES	1	0	0
OXA-23	2	0	1
OXA-58	2	0	0
OXA-24/40	0	0	1
OXA-51	0	0	0
NEG	12	0	10

Some strains/samples had several target genes.

Results

- Results: see Table 2
- High sensitivity: 100% with isolates and samples
- Specificity: 100% with strains
- Lower specificity (60%) with clinical samples (directly without culturing) can be explained by the higher sensitivity of these methods compared with the culture:
 - 6/7 of false positive samples were from patients who had matching strains and genes in other samples
 - The only "real" false positive was OXA-40 which was only detected with Eazyplex SuperBug Complete B, and the matching strain/gene could not be identified in other samples of the same patient

Table 2. Results

		eazyplex® SuperBug complete B	NOVODIAG CarbaR+	AMPLIDIAG CarbaR+MCR
Bacterial isolates	Sensitivity	100%	100%	100%
	Specificity	100%	100%	100%
	True POS	21	28	27
	True NEG	18	13	12
Clinical samples	False POS	0	0	0
	False NEG	0	0	0
	Sensitivity	100%	100%	100%
	Specificity	59%	63%	63%
	True POS	25	25	25
	True NEG	10	10	10
	False POS	7	6	6
	False NEG	0	0	0

Conclusions

- Sensitivity: 100% (isolates and samples)
- Specificity: 100% (isolates) ca. 60% (samples): this can be explained by the superiority of molecular methods to our reference culture method
- All three commercial CPE tests were easy-to-use and reliable methods for the detection of carbapenemase genes from bacterial colonies and clinical samples

Acknowledgements

We thank Mira Id, Milla Kulin, and Janita Toikkanen for technical assistance.

Disclosures

The authors have no conflict of interest.

References

Pasanen T, Koskela S, Mero S, Tarkka E, Tissari P, Vaara M, Kirveskari J. PLoS One. 2014;9(1):e85854. doi: 10.1371/journal.pone.0085854. Rapid molecular characterization of *Acinetobacter baumannii* clones with rep-PCR and evaluation of carbapenemase genes by new multiplex PCR in Hospital District of Helsinki and Uusimaa.